

BC 760 & BC 780 Technical Specifications



Principles

WBC (MG/Neu/Mon/Lym/Eos/Bas), NRBC/RET*, PLT-H/PLT-O*/IPF, SF Cube*
Cell Analysis Technology *S: Scatter; F: Fluorescence; Cube: 3D analysis

RBC, PLT

Focusing Flow-DC Impedance Method

HGB

Colorimetric method

ESR

Photometric method

Number of measuring parameters (whole blood): 109

Number of reportable parameters: 41*

WBC Bas#, Bas%, Neu#, Neu%, Eos#, Eos%, Lym#, Lym%, Mon#, Mon%, IMG#, IMG%, RET%, RET#, RHE, IRF*, LFR*, MFR*, HFR*, RBC HGB, MCV, MCH, MCHC, RDW-CV, RDW-SD, HCT, NRBC#, NRBC%, PLT, PLT-H, PLT-O*, MPV, PDW, PCT, P-LCR, P-LCC, IPF, ESR*

Number of research parameters: 68*

Number of measuring parameters (body fluid): 18

Number of reportable parameters: 7

WBC-BF, TC-BF#, MNe, MN%, PMNe, PMN%, RBC-BF

Number of research parameters: 11

Data storage capacity

Up to 150,000 results including numeric and graphical information*

Analysis Mode

Sample Type	Analysis Mode
Whole blood	CBC, CBC + DIFF, CBC + DIFF+RET*, CD + ESR, CDR + ESR*, CD/WBC-3X, CDR/PLT-5X*, and other modes
Predilute	CBC, CBC + DIFF, CDR*, and other modes
Body fluid	CBC + DIFF

Physical Specifications

*Dimensions: 840D x 655W x 600H mm
Weight: ≤70.6Kg

Voltage: 100V-240V~ (±10%)

Frequency: 50Hz/60Hz (±1Hz)

Power input: 600VA

External output

LANx1, USB x 4 (Specifications: DC 5V; 500mA; USB2.0 x 3; USB3.0 x 1)

Normal Operating Environment

Ambient temperature: 10 °C~ 35 °C

Relative Humidity:

Atmospheric pressure:

70.0kPa ~ 106.0kPa*

*Note: Required altitude for normal operation:
-400m ~ +3000m*

Sample volume

CD (whole blood): 25ul

CD+ESR (whole blood): 160ul

Predilute: 20ul

Throughput

CD 80t/h CDR 45t/h CD+ESR 40t/h

Performance

Parameter	Linearity Range	Precision	Carryover
WBC	0-500×10 ⁹ /L	≤2.5% (≥4.51×10 ⁹ /L)	≤1.0%
RBC	0-8.60×10 ¹² /L	≤1.5% (≥3.5×10 ¹² /L)	≤1.0%
HGB	0-260g/L	≤1.0% (110-180g/L)	≤1.0%
HCT	0-75%	≤1.5% (30%-50%)	≤1.0%
PLT*	0-5000×10 ⁹ /L	≤1.5(SD) (≤20×10 ⁹ /L) [^] ≤2.5% (≥100×10 ⁹ /L) [^]	≤1.0%
RET*	0-0.8×10 ¹² /L	≤15% (RBC ≥ 3.00×10 ¹² /L) RET%: 1.00% ~ 4.00%	≤1.0%
ESR		≤1.8(SD)(0-20mm/h)	≤1.0%

Note: Applicable only to CDR/PLT-O 5x and CR/PLT-O 5x models

Items marked with an asterisk (*) apply only to BC-780

BC 760 & BC 780 Advanced Automated Hematology Analyzer

ESR Testing - Above & Beyond



VWR.COM

Avantor Performance Materials India Limited

17th Floor, Building No. 5, Tower C, DLF Cyber City, Phase - III, Gurgaon - 122002, Haryana
Tel: +91 124 4656700 email: csdiagnostics@avantorsciences.com

Marketed by:

Abaxis Healthcare India Pvt. Ltd.

Peringavu, Thrissur, Kerala-680008, Phone: +91 487 2323537


mail: info@abaxishealthcare.com www.abaxishealthcare.com



abaxis healthcare

BC 760 & BC 780

Automated Hematology Analyzer





Above your expectations

SF Cube fluorescent technology allows reliable counting and differentiation of abnormal samples.

More refined and reliable cell differentiation

3D fluorescent analysis technology allows reliable differentiation of immature and other abnormal cells, such as immature granulocytes (IMGs), reticulocytes (RETs*), and immature platelet fraction (IPF)



More reliable measurements for low-value samples

The BC-760 & BC-780 3D fluorescence analysis platform is designed with multiple counting WBC-3x and PLT-O 5x analysis modes to help ensure higher reliability for low-value WBC and PLT samples. In addition, the PLT de-aggregation function can reduce the cumbersome review work.

More comprehensive alarm messages for abnormalities

The analyzer provides a detailed list of over 40 prompt messages, including WBC message, RBC message, and PLT message. This allows laboratory technicians to intuitively and quickly identify abnormal samples and proceed further with the samples in a timely manner. This in turn helps to avoid missed diagnosis of blood disease and false reports

Parameter	Result	Unit	Parameter	Result	Unit
WBC	15.54	10 ⁹ /L	RBC	4.99	10 ¹² /L
Hgb	15.30	g/L	HCT	5.18	%
Hct	0.35	%	MCV	102.8	fL
MCH	0.41	pg	MCHC	22.8	g/dL
MCHC	0.29	g/dL	MCN	21.7	10 ⁹ /L
RDW	0.02	%	RDW-CV	1.30	%
RDW-CV	0.41	%	RDW-SD	17.7	fL
RDW-SD	0.02	%	RET	0.004	%
RET	0.004	%	RET-H	0.01	%
RET-H	0.001	%	IPF	0.0	%
IPF	0.001	%	PLT	14.8	10 ⁹ /L
PLT	14.8	10 ⁹ /L	PLT-H	0.00	%
PLT-H	0.00	%	ESR	0.0	mm/h
ESR	0.0	mm/h			

WBC Message


PLT Message

ESR Message

Beyond your expectations


Limitations of traditional PLT counting

In the traditional impedance method, PLTs are subject to interferences that may lead to falsely high or falsely low results (as shown in the figure). Once an error report is generated, it will directly affect the judgment and decision-making of clinicians. The results reported at the clinical decision level are related to patient safety. Therefore, accurate PLT results are critical in clinical practice.



Optical PLT-H in every CD test

In order to solve the above problem, we have developed a brand new parameter PLT-H. It combines small PLTs from the conventional impedance method and large PLTs from the optical method. The solution can resist the interferences in conventional PLT detection without requiring extra reagents.



CD + ESR in one test provide reliable ESR results with greater ease

The BC-700 series integrates an automatic ESR module in a hematology analyzer. It can also generate both CBC & ESR results in one test within 1.5 min.

In addition, it saves the costs that would otherwise be incurred for the purchase, maintenance, consumables, and storage space of a separate ESR analyzer. Compared with the traditional Westergren method, this method performs better in quality traceability, repeatability, speed, safety, and level of automation.

